

AW-KS-Set-ICT: Extraction Tool for Receptacles of series GKS-040 / 050 / 075 / 100

Using this Extraction Tool wrongly set and damaged Receptacles can be easily removed within seconds.

The Tool cuts itself into the outer wall of the Receptacle when screwed in. The spring-loaded weight enables extraction of the Receptacle out of the assembly hole without breakage or further damage.


Combining Extraction Tool, Adapter and Insert accordingly, Receptacles of the series GKS-040/GKS-050/GKS-075 and GKS-100 can be removed from the Probe Plate.

The handy Tool-Box AW-KS-ICT offers all components which are needed for simple removing:

- 1 x Extraction Tool (pre-assembled with Standard Spring)
- 1 x Adapter for Bits for extraction of KS-100
- 1 x Adapter for Bits for extraction of KS-075, KS-050 und KS-040
- 3 x Bit for KS-100
- 3 x Bit for KS-075
- 3 x Bit for KS-050
- 3 x Bit for KS-040
- 1 x Spring with higher spring force
- 1 x Hexagon Key SW0.9 for Threaded Pin

Order Description:

AW-KS-Set-ICT



AW-KS-Set-ICT

Please contact us for further information.
Prices and delivery time on request.
Technical changes possible without prior notification.